

User's Manual and Instructions

Express Cloning Checker Kits

Catalog Number:

K501 1200	K5012100	K5013200
K501 1200-1	K5012100-1	
K501 1200-2	K5012100-2	

Introduction

Express Cloning Checker kit enables you to perform large-scale screening of up to hundreds of colonies, in half a day. In addition, it can be used to check insert sizes in constructed cDNA libraries without the preparation of plasmid DNA. In the same way, plasmid DNA can be quickly analyzed by directly using liquid bacterial culture, glycerol strain stock without the extraction of plasmid DNA. This is especially useful in the case of insufficient ligation generating very high background, or the cloning of very large inserts.

Description

Express Cloning Checker Kits provide two highly efficient methods of identifying recombinant colonies after transformation without time-consuming plasmid preparation. By using the large-scale screen method in Kit I, all you need to do is to transfer partial colonies of E. coli bacteria directly from the transformation plates into the solutions provided by the kits and immediately run an agarose gel electrophoresis. Within one to two hours, you will know exactly which are recombinant colonies, instead of spending two days using time and labor-intensive analysis methods such as plasmid DNA miniprep or colony hybridization. In case exact size or orientation of inserts is of special concern, there is a second choice. The enzyme digestion method of Kit II, in which a unique solution is used to release the plasmid DNA from host bacteria cells without affecting following reactions of restriction enzymes. So, one can perform a regular restriction enzyme analysis in a single tube by directly using bacteria without preparing plasmid DNA. With the combination of these two methods found in Kit III, one can use the first method to perform quick and large-scale screening in the expressway and use the second method to confirm the recombinants.

Features

- Rapid and efficient identification of recombinant colonies within 1 hour
- Large-scale screening up to hundreds of colonies in short time
- Quick check of insert sizes in constructed cDNA libraries
- Restriction enzyme analyzing of recombinants directly from bacteria without preparation of plasmid DNA
- Compatible with most commonly used E. coli. bacteria from plate colonies, liquid culture or glycerol stock
- Replacing time-consuming DNA miniprep and restriction enzyme digestions or colony PCR

Kit contents and storage

There are three types of kit. Check which kit you have.

Kit I Cat # K5011200 (200 reactions) Previous Cat# 035301 Shipment condition: Room Temperature

Contents	Volume	Storage
Red solution	1.0 ml	Room temperature
Yellow solution	1.0 ml	Room temperature
Supercoiled DNA marker* (size: 3.0, 3.9, 5.4, 7.0, 10 kb)	2 µg	Room temperature
Screening plates	1	Room temperature

Kit II Cat # K5012100 (100 reactions) Previous Cat# 035302 Shipment condition: Room Temperature

Contents	Volume	Storage
Green solution	0.8 ml	Room temperature
Blue solution	0.2 ml	Room temperature

Kit III Cat # K5013200 (200 reactions) Previous Cat# 035303 Shipment condition: Room Temperature

Contents	Volume	Storage
Red solution	0.5 ml	Room temperature
Yellow solution	0.5 ml	Room temperature
Green solution	0.8 ml	Room temperature
Blue solution	0.2 ml	Room temperature
Supercoiled DNA marker* (size: 3.0, 3.9, 5.4, 7.0, 10 kb)	2 µg	Room temperature
Screening plates	1	Room temperature

- Supercoiled DNA marker has been freeze-dried for room temperature storage, add 20 µl sterile water and mix by vortex and spinning briefly several rounds before use. The resulted marker concentration is 100 ng/µl and should be stored at -20 °C.

Brief protocol:

A. Large scale screen method (Kit I and Kit III):

1. Pick up bacteria colonies and stir in 5 µl of Red solution.
2. Add 5 µl of Yellow solution.
3. Run gel and check results.

B. Enzyme digest method (Kit II and Kit III):

1. Pick up bacteria colonies and stir in 8 µl of Green solution.
2. Heat at 100° C for 30 seconds.
3. Add 2 µl of enzyme mix and incubate at appropriate temperature for 30 minutes.
4. Add 2 µl of Blue solution.
5. Run gel and check results.

Special notice for using the kit:

- Use freshly prepared culture plates
- Spread less 300 bacterial colonies per plate
- Post-stain gel (no EB in gel and running buffer)
- Add 20 µl sterile water in Supercoiled DNA marker tube and mix it before use. Store at -20 °C.

Detailed protocol

A. Supercoiled Plasmid Analysis (for Kit I and Kit III)

After overnight culture (16- 24 hours) of transformation reaction:

1. Pour a 0.8 to 1.2 % agarose gel with thin comb (thickness less than 1 mm recommended). Do not submerge the gel in running buffer before sample loading.
2. Mark 10 to 30 well separated colonies with numbers on the bottoms of transformation plates¹.
3. Add 5 µl of Red solution per well of screening plate supplied by the kit or micro-tubes. The screen plate can be used repeatedly after wash.
4. Pick up partial bacterial colony² (50 % to 25% if a colony size around 1.5 - 2.5 mm in diameter) of each selected colony with pipette tips. Stir attached bacteria 2 or 3 times into Red solution in wells or tubes from step 2. If the whole colony has been picked up, place the plate back to 37° C and incubate for a desired period after this step is done.
5. Prepare DNA marker. a) Use supercoiled DNA marker provided in kit. Add 20 µl of sterile water into DNA marker and mix by vortex and spinning briefly several rounds before the first use. Take 2 µl of the marker and add it into 5 µl of Red

¹ Use freshly prepared plates to culture bacterial colonies. To get well-separated and well-grown colonies, spreading of less 300 colonies per plate is recommended. Old bacterial colonies in plates, which are stored or cultured for several days, may produce weaker DNA bands.

² Take 1 to 3 µl of bacterial culture or plasmid stock strain instead of colonies if bacterial culture of a clone is available.

solution. b). Use your own marker. Take 10 - 30 ng (in 2 μ l volume) of vector DNA without insert and add it into 5 μ l the Red solution.

6. Add 5 μ l of Yellow solution to each well or tube containing bacteria or DNA marker.
7. Mix samples by pipeting up and down 2 or 3 times and load them directly into loading wells of agarose gel³ before setting the gel in running buffer.

Option: if using micro-tubes, but not supplied screening plate from step 3, one may load the samples in submerged gel in running buffer as an usual electrophoresis by following steps: 1) briefly mix samples after adding of Yellow Solutions; 2) add 4 μ l of Blue buffer to each reaction; 3) vortex samples for about 15 seconds before loading samples. These extra steps can increase target DNA signals and avoid defused bands generated from a none-fresh gel. So, it is only recommended to use if checking low-copy number plasmid clones or screening large quantity of colonies (more than 40).

8. Set the gel with loaded samples in running buffer (1x TAE or TBE) and begin electrophoresis. Make sure running buffer gets into every loading well without bubbles. After the red dye band⁴ reaching 3-4 cm from loading well (cloning of small inserts less than 300 bp needs a little longer running), stop running and visualize the DNA by staining the gel in 0.2 - 0.5 μ g/ml ethidium bromide⁵.
9. Examine the gel by UV transillumination. Due to low-copy number of plasmid DNA in some bacterial strain, photography of gel may need to use high aperture (4.5) and long exposure (2 to 5 seconds). Recombinant colonies will migrate slower in the gel than the vector without insert⁶. Now go back the original plates and select corresponding recombinant colonies to culture or stock for further analysis.

B. Restriction Enzyme Digestion Analysis (for Kit II and Kit III to check size and orientation of inserts or restriction site information)

1. Mark 10 to 30 well separated colonies with numbers on the bottoms of transformation plates¹.
2. Add 8 μ l of green solution in each microcentrifuge tube or wells of 96-well plate.
3. Pick up partial bacterial colony² (50% to 25% if a colony size around 1.5 - 2.5 mm in diameter) of each selected colony with small pipet tips⁷ and stir the attached bacteria 2 or 3 times into Green solution in tubes or wells from step 2. If the whole colony has been picked up, place the plate back to 37° C and incubate for a desired while after the step is done.
4. Heat the samples at 100 ° C for 30 seconds in boiled water or by using PCR machine.
5. After the samples cool down to room temperature, add 2 μ l of digest mix containing 1 μ l of 10 x digest buffer and appropriate restriction enzymes (1 - 5 unit per reaction).
6. Incubate the samples at appropriate temperature for 10 - 30 min⁸.
7. Pour a 0.8 to 1.5 % agarose gel³ with thin comb (thickness less than 1 mm recommended).
8. Add 2 μ l of Blue solution in each reaction and mix them.
9. Load samples in wells of the agarose gel. Finally, load linearized DNA size marker (not supplied in the kit) for DNA mobility comparison.
10. Begin electrophoresis in running buffer (1x TAE or TBE). Make sure running buffer gets into every loading well without bubbles. After the first dye (blue band) reaching about 3 - 4 cm length from loading well of the gel, stop running and visualize the DNA by staining the gel in 0.2 - 0.5 μ g/ml ethidium bromide⁴.
11. Examine the gel by UV transillumination. Due to low copy number of plasmid DNA in some bacterial strain, photography of gel may need to use high aperture (4.5) and long exposure (2 - 5 seconds). Now go back to the original plates and select corresponding recombinant colonies to culture or stock for further analysis.

³ It is recommended to use freshly prepared gel (within 30 minutes after pouring off). Over-dried gel will generate diffused bands. Use 1.5% agarose gel if checking small inserts (less than 300 bp).

⁴ Red dye migrates through agarose gel slower than bromophenol blue at approximately the same rate as linear double-stranded DNA 1.2 kb in length.

⁵ It is strongly recommended using post-staining of gel and fresh staining buffer, which is made from EB stock solution. Adding EB in electrophoresis buffer or in gel will result in inconsistent migration of plasmid DNA.

⁶ There may be four bands observable (like lane 5, Fig 1) for each colony if picking up large piece of colonies and gel running time long enough. The fastest moving DNA bands, which are bigger than 2.5 kb, are strong and informative. Sometimes non-nucleic acid material and RNA will show up like diffused bands on bottom of gel, so comparison should be made between the informative bands of colony plasmid and that of control markers.

⁷ Avoid using of toothpicks to pick up colonies at this step because it absorbs the green solution.

⁸ Add more enzymes will shorten the incubation period. Over-incubation may result in digesting of genomic DNA and increasing background.

Analysis of recombinant colonies with 1% agarose gel

MC 1 2 3 4 5 6




Fig 1. Large scale screening analysis of recombinant colonies with different size of inserts which were demonstrated in lane 1 through 6 with comparison of the control vector without insert DNA (lane C, 5.4 kb). Lane M is *Hind* III/ lambda DNA marker.

M 1 2 3 4 5 6



Fig 2. Enzyme digest analysis of the extracted plasmid DNA of these recombinant colonies in Fig 2 after minipreparation showing that they contains 0.9, 1.2, 1.7, 1.9, 2.2 and 2.5 kb insert respectively. At least 200 bp insert can be distinguished by the large scale screen method.

1 2 3 4 5 6 7 8 9 10 M




Fig 3. Large scale screening of recombinant colonies. Partial colonies were picked up directly from overnight cultured transformation plates and treated with red and yellow solutions before electrophoresis. Three recombinants with 2.5 kb insert (lanes 3, 5 and 8) are identified. Lane M is the supercoiled DNA marker supplied in the kit.




Fig 4. Enzyme digest analysis of recombinant colonies. Recombinant colonies were picked up directly from overnight cultured transformation plates in green solution and followed by restriction enzyme digestion (*Hind* III/*Xba* I, lanes 1 and 2; *EcoR* I lanes 3 and 4) for 30 minutes at 37 °C. Lane M is *Hind* III/ lambda DNA marker.