



Sample Buffer without detergents

Antibody dilution buffer for immunohistochemistry

Product description

pH-value: pH 7.2 ± 0.2

Preservative: contains 0.1% ProClin® 300

Storage: shelf life at –20 °C, 1 year – tolerates repeated freezing and thawing cycles
shelf life at 2 – 8 °C, 6 months

For research use only, not for diagnostic use

Available

Package sizes:	50 ml	order number 106 050
	125 ml	order number 106 125
	500 ml	order number 106 500

Instructions for use

Sample Buffer is a dilution buffer for specimen and antibodies for use in immunoassays. Antibodies, coupled to alkaline phosphatase or peroxidase, can be diluted and used in *Sample Buffer* without problems, because *Sample Buffer* is free of phosphate.

Sample Buffer without detergents contains neither Tween nor other detergents. Thus it is for use especially in **immunohistochemistry**, where detergents sometimes can make problems in detection.

Immediately before using the buffer should be mixed thoroughly.

Sample Buffer is ready-to-use. Repeated freezing and thawing is possible.

Primary and secondary antibodies can be diluted directly in *Sample Buffer*.

Sample Buffer is designed for trouble free assays without interference. If you have background, we recommend using *LowCross®-Buffer* rather than *Sample Buffer*.

For further information please visit www.candor-bioscience.com.

ProClin is a registered trade mark of Rohm and Haas Company.
LowCross-Buffer is a registered trade mark of CANDOR Bioscience GmbH.