


Wnt/ β -catenin Signaling Research

TCF/LEF Dual Reporter Vector (Cat. #9038)

Precisely **monitor Wnt pathway activation** with ScienCell's TCF/LEF Dual Reporter Vector, featuring a **Firefly luciferase** under TCF/LEF control and a **Renilla luciferase** for normalization. Includes a **negative control vector** for accurate background correction. Suitable for stem cell research, oncology, developmental biology, and drug discovery.



Related Wnt Protein Products

Product Name	Product Code	Catalog #
Recombinant Human Wnt-3a Protein, E. coli Derived	rhWnt3AE	109-03AE
Recombinant Human Wnt-3a Protein Fragment 174-240	rhWnt3AEF10	109-03AEF10
Recombinant Human Wnt-3a Protein Fragment 198-223	rhWnt3AEF11	109-03AEF11
Recombinant Human Wnt-3a Protein Fragment 19-87, 294-352	rhWnt3AEF12	109-03AEF12
Recombinant Human Wnt-3a Protein Fragment 19-87	rhWnt3AEF2	109-03AEF2
Recombinant Human Wnt-3a Protein Fragment 216-352	rhWnt3AEF4	109-03AEF4
Recombinant Human Wnt-3a Protein Fragment 19-105, 241-352	rhWnt3AEF5	109-03AEF5
Recombinant Human Wnt-3a Protein Fragment 223-352	rhWnt3AEF6	109-03AEF6
Recombinant Human Wnt-3a Protein Fragment 223-269	rhWnt3AEF7	109-03AEF7




Discover ScienCell's Cell-Based Assays

UniFectagen (Cat. #0983)

Enables **high-efficiency transfection** of plasmid DNA, siRNA, and viral vectors into mammalian cells with **minimal cytotoxicity**. It delivers **strong gene expression** across cell lines such as HEK293T and CHO, in a ready-to-use, serum-compatible format for seamless workflow integration.

Luciferase Assay Kit (Cat. #9048)

Our Dual Luciferase Assay Kit enables easy measurement of Firefly and Renilla luciferase in the same sample with high precision. Separate substrates ensure clean, accurate results every time—ideal for **reporter gene studies**, **promoter analysis**, and **signaling pathway research**.



Firefly Assay Enhancer (Cat. #7208)

ScienCell's Firefly Assay Enhancer **boosts luminescence intensity** and **signal duration** for stronger, more consistent Firefly Luciferase Assay results. Compatible with Luciferase Assay Kit (Cat. # 9048) and available in a **ready-to-use** format for easy workflow integration.

Human Cancer Promoter Methylation qPCR Assay Kit (Cat. #9058)

Enables **sensitive and accurate detection** of promoter methylation in key **cancer-associated genes (RASSF1A, CDKN2A/p16, and GSTP1)** using high-specificity primers and a built-in 1:1 methylated DNA control for precise quantification, supporting applications in **cancer biomarker discovery**, **epigenetic profiling**, and **translational research**.

