

Specialty Chemicals Decontaminants & Cleaners

Decontaminants & Cleaners

AbSolve Glassware Cleaner

- Easily and inexpensively removes serious R or DNA contamination (autoclaving does not destroy DNA). Simply soak plastic and glassware in a 2% solution for 30 minutes, rinse and use.
- Ideal for prewashing plastic tubes for PCR. AbSolve™ prevents DNA contamination of storage and reaction vessels from exogeneous sources.
- Use to clean glassware for electrophoresis. Glass plates rinse easily and AbSolve leaves no residue to interfere with gel polymerization or silver staining.
- Safe, non-abrasive AbSolve does not etch glass, contain strong acids, or emit toxic fumes.

COUNT-OFF Liquid Concentrate

- Ideal, all-purpose decontaminant: safe, efficient and economical.
- Quickly and effectively cleans up even the most persistent leftovers: proven most effective cleanser for removing radioactive residues (³H, ¹⁴C, ³²P and ¹²⁵I), stopcock and vacuum greases, lanolin and petroleum jelly, dried blood and serum, fatty and amino acids, protein complexes, and polymer films and other stubborn residues.
- Stable under extreme temperature fluctuations (-50 °C to 150 °C).
- Safer than strongly acidic cleaners, such as chromic acid: solutions will not produce toxic gases from substrates containing ¹⁴C, ¹³¹I, ³⁵S, or ³⁶Cl.
- Easy to use: adjust concentration, soaking time and temperature as needed to speed up decontamination or for particularly stubborn substances such as dried blood.

COUNT-OFF Surface Cleaner (Trigger Spray)

- Quickly and safely decontaminates small radioactive spills from benches, shields, and appliances: also efficiently cleans instrument housings, hood corners, centrifuge cups and heads, and LSC counter mechanisms.
- Traps and suspends radioactive particles in foam: makes clean-up easy and reduces the likelihood of spreading contamination to hands and clothing.
- Removes both ionic and non-polar radioactivity.
- Cleans-up stubborn grease, resins, blood and wax, even from rough surfaces.

Pico-Kleen N

Pico-Kleen N is a concentrated liquid detergent developed as an all-purpose cleaning agent and radioactive decontaminant. Pico-Kleen N is essentially neutral, in which mildness is combined with an outstanding surfactant formulation resulting in an effective cleaning action. Depending on the type of radioactive contaminated surface and isotopes, Pico-Kleen N will effectively remove the contamination when properly applied.

Product	Product No.	Size
AbSolve™ Glassware Cleaner	6NE9711	1 x 1 L
COUNT-OFF™ Liquid Concentrate	6NE9422	1 x 2.5 L
	6NE9427	4 x 2.5 L
COUNT-OFF Surface Cleaner (Trigger Spray)	6NE942T	6 x 22 oz. (650 mL) pump bottles
Pico-Kleen N	6013819	2 x 5 L