

Data Sheet

Product Name: Frozen Tissue Array**Catalog No.:** T6234700-5**Lot No.:** B907113**Species:** ☒ Human ☐ Mouse ☐ Rat ☐ Monkey (Rh) ☐ Guinea Pig ☐ Porcine
☐ Bovine ☐ Hamster ☐ Dog ☐ Monkey (Cy) ☐ Rabbit ☐ Plant**Tissue Type:** ☒ Normal ☒ Adult ☐ Fetal ☐ Tumor ☐ Disease ☐ Cell line**Tissue Array Diagram:**

	1	2	3	4	5	6	7	8	
A	☐	☐	☐	☐	☐	☐	☐	☐	
B	☐	☐	☐	☐	☐	☐	☐	☐	
C	☐	☐	☐	☐	☐	☐	☐	☐	
D	☐	☐	☐	☐	☐	☐	☐	☐	
E	☐	☐	☐	☐	☐	☐	☐	☐	

Components:

- ☐ 2 Tissue Array Slides ☒ 5 Tissue Array Slides per package
- Certificate of Analysis

FOR IN VITRO RESEARCH USE ONLY

APPROVED BY: _____



Data Sheet

Product Name: Frozen Tissue Array**Catalog No.:** T6234700-5**Lot No.:** B907113**Donor Information**

Position	Tissue Type	Age	Gender	Diagnosis
A1, A2	Adrenal	47	M	Normal
A3, A4	Ovary	20	F	Normal
A5, A6	Pancreas	44	M	Normal
A7, A8	Thyroid	25	M	Normal
B1, B2	Brain, Cerebrum	70	F	Normal
B3, B4	Brain, Cerebellum	49	M	Normal
B5, B6	Lung	83	F	Normal
B7, B8	Spleen	30	M	Normal
C1, C2	Uterus	41	F	Normal
C3, C4	Cervix	30	F	Normal
C5, C6	Breast	64	F	Normal
C7, C8	Placenta	24	F	Normal
D1, D2	Heart	87	F	Normal
D3, D4	Skin	78	F	Normal
D5, D6	Skeletal Muscle	79	M	Normal
D7, D8	Kidney	50	M	Normal
E1, E2	Stomach	56	M	Normal
E3, E4	Small Intestine	43	F	Normal
E5, E6	Liver	30	M	Normal
E7, E8	Salivary Gland	87	F	Normal

FOR IN VITRO RESEARCH USE ONLY**APPROVED BY:** 