

Data Sheet

Product Name: Frozen Tumor and Normal Tissue Array

Catalog No.: T6235700-1

Lot No.: B712100

Species: ☒ Human ☐ Mouse ☐ Rat ☐ Monkey (Rh) ☐ Guinea Pig ☐ Porcine
☐ Bovine ☐ Hamster ☐ Dog ☐ Monkey (Cy) ☐ Rabbit ☐ Plant

Tissue Type: ☒ Normal ☐ Adult ☐ Fetal ☒ Tumor ☐ Disease ☐ Cell line

Tissue Array Diagram:

	1	2	3	4	5	6	7	
A	☐	☐	☐	☐	☐	☐	☐	
B	☐	☐	☐	☐	☐	☐	☐	
C	☐	☐	☐	☐	☐	☐	☐	
D	☐	☐	☐	☐	☐	☐	☐	

Components:

1. 1 Tissue Array Slides per package
2. Certificate of Analysis

FOR IN VITRO RESEARCH USE ONLY

APPROVED BY: _____



Data Sheet

Product Name: Frozen Tumor and Normal Tissue Array

Catalog No.: T6235700-1

Lot No.: B712100

Donor Information:

Position	Tissue Type	Age	Sex	Pathological Diagnosis	Differentiation	TNM or Stage
A1	Brain	70	F	Normal		
B1	Brain Tumor	36	F	Oligodendroglioma	N/A	Stage III
A2	Breast	40	F	Normal		
B2	Breast Tumor	47	F	Invasive ductal carcinoma	N/A	Tumor size unknownNoM0
A3	Colon	87	F	Normal		
B3	Colon Tumor	70	M	Adenocarcinoma, mucuous	Moderately	T2NoM0
A4	Skeletal Muscle	79	M	Normal		
B4	Skeletal Muscle Tumor	50	M	Rhabdomyosarcoma	Poorly	T3NoM0
A5	Kidney	44	M	Normal		
B5	Kidney Tumor	37	M	Renal Cell Carcinoma	Moderately	T3NoM1
A6	Liver	64	M	Normal		
B6	Liver Tumor	44	M	Hepatocellular Carcinoma	N/A	T3NoM0
A7	Lung	83	F	Normal		
B7	Lung Tumor	70	M	Adenocarcinoma	Moderately	Tumor size unknownNoM0
C1	Pancreas	86	F	Normal		
D1	Pancreas Tumor	53	M	Adenocarcinoma	Poorly	Tumor size unknownNoM0
C2	Prostate	50	M	Normal		
D2	Prostate Tumor	66	M	Adenocarcinoma	N/A	Gleason 4+3=7
C3	Skin	61	F	Normal		
D3	Skin Tumor	48	M	Cacinoma, sweat gland	N/A	T1NoM0
C4	Small Intestine	70	F	Normal		
D4	Small Intestine Tumor	68	M	Malignant Mesenchymoma	Well	T2NoM1
C5	Stomach	56	M	Normal		
D5	Stomach Tumor	54	M	Adenocarcinoma, ulcer	Moderately	T2NoM0
C6	Ovary	37	F	Normal		
D6	Ovary Tumor	54	F	Cystadenocarcinoma, Serous	Poorly	T2NoM0
C7	Uterus	68	F	Normal		
D7	Uterus Tumor	55	F	Adenocarcinoma	Poorly	Tumor size unknownNoM0

FOR IN VITRO RESEARCH USE ONLY

APPROVED BY: _____

