


細胞影像觀察攝影 傷口癒合分析 最佳實驗工具



指定商品 **第2件半價** 優惠


優惠時間：2019.Sep.2 - Oct.16




超薄底設計，符合 No. 1.5 蓋玻片厚度，影像清晰不變形！



細胞培養 → 免疫染色 → 影像觀察，一氣呵成！



ibiTreat 表面處理，適合各類貼附型細胞生長！



Culture-Insert 提供大小一致的細胞間隙，實驗結果更精確

應用實驗

細胞影像觀察



- ① μ-Dish 35 mm
- ② μ-Slide 2 well
- ③ μ-Slide 4 well
- ④ μ-Slide 8 well
- ⑤ 3, 8, 12 well Chamber
- ⑥ μ-Slide 2 x 9 well


優惠產品

包裝

貨號

μ-Dish 35 mm, low ibiTreat	60	IB-80136
μ-Dish 35 mm, low ibiTreat Grid-500	60	IB-80156
μ-Slide 2 well ibiTreat	15	IB-80286
μ-Slide 2 well Ph+ ibiTreat	15	IB-80296
μ-Slide 4 well ibiTreat	15	IB-80426
μ-Slide 4 well Ph+ ibiTreat	15	IB-80446
μ-Slide 8 well ibiTreat	15	IB-80826
μ-Dish 35 mm, high Glass Bottom Grid-50	30	IB-81148
μ-Dish 35 mm, high ibiTreat	60	IB-81156
μ-Dish 35 mm, high Glass Bottom	60	IB-81158
μ-Dish 35 mm, high ibiTreat Grid-500	60	IB-81166
μ-Dish 35 mm, high Glass Bottom Grid-500	30	IB-81168
3 well Chamber, removable	15	IB-80381
8 well Chamber, removable	15	IB-80841
12 well Chamber, removable	15	IB-81201
μ-Slide 2 x 9 well ibiTreat	15	IB-81806
25 Culture-Inserts for self-insertion	25	IB-80209
Culture-Insert 24 ibiTreat	3	IB-80241
Culture-Insert in μ-Dish 35 mm, high ibiTreat	30	IB-81176
Culture-Insert 3 well in μ-Dish 35 mm, high ibiTreat	30	IB-80366
25 Culture-Inserts 3 well for self-insertion	25	IB-80369
Culture-Insert 4 well in μ-Dish 35 mm, high ibiTreat	30	IB-80466
25 Culture-Inserts 4 well for self-insertion	25	IB-80469

傷口癒合/細胞遷移試驗



- ① Culture-Insert
- ② Culture-Insert 24
- ③ Culture-Insert 3 well
- ④ Culture-Insert 4 well

