



F-actin 螢光染色試劑

精確掌握細胞內 F-actin 分布情形與表現量變化







CytoPainter Phalloidin F-actin Staining Kits and Reagents



Phalloidin 為一種帶有毒性的胜肽，具有與纖維形肌動蛋白 (filamentous actin, F-actin) 專一性結合的特性。Abcam 「CytoPainter Phalloidin」 為標定有螢光染料的 Phalloidin 試劑產品，能夠協助研究人員精確掌握細胞內 F-actin 的分布情形與表現量變化。

- 適用於多種樣本類型 — Fixed and permeabilized tissue sections、Adherent/Suspension cells 等。
- 螢光標定選擇多樣化 — FITC、Rhodamine、iFluor 405、iFluor 488、iFluor 555、iFluor 647 等。
- 兩種產品包裝型式 — Kit 內含 10 倍濃縮螢光染劑及稀釋液；Reagent 僅含 10 倍濃縮螢光染劑。

Product Name	Ex/Em	Pack Size	Cat. No.
 F-actin Staining Kit - Blue Fluorescence - Cytopainter	350/450	500 tests	ab112124
 F-actin Staining Kit - Green Fluorescence - Cytopainter	490/525	500 tests	ab112125
 F-actin Staining Kit - Orange Fluorescence - Cytopainter	550/575	500 tests	ab112126
 F-actin Staining Kit - Red Fluorescence - Cytopainter	594/610	500 tests	ab112127
 Phalloidin-iFluor 350 Reagent	353/442	300 tests	ab176751
 Phalloidin-iFluor 405 Reagent	400/421	300 tests	ab176752
 Phalloidin FITC Reagent	492/518	300 tests	ab235137
 Phalloidin-iFluor 488 Reagent	493/517	300 tests	ab176753
 Phalloidin-iFluor 514 Reagent	520/547	300 tests	ab176754
 Phalloidin-iFluor 532 Reagent	542/558	300 tests	ab176755
 Rhodamine Phalloidin Reagent	546/575	300 tests	ab235138
 Phalloidin-iFluor 555 Reagent	556/574	300 tests	ab176756
 Phalloidin-iFluor 594 Reagent	590/618	300 tests	ab176757
 Phalloidin-iFluor 633 Reagent	634/649	300 tests	ab176758
 Phalloidin-iFluor 647 Reagent	650/665	300 tests	ab176759
 Phalloidin-iFluor 680 Reagent	683/700	300 tests	ab176760
 Phalloidin-iFluor 700 Reagent	692/708	300 tests	ab176761
 Phalloidin-iFluor 750 Reagent	752/778	300 tests	ab176762
 Phalloidin-iFluor 790 Reagent	787/808	300 tests	ab176763



Phalloidin staining protocol

