

細胞螢光染劑

- ✓ 非抗體螢光探針，可專一性追蹤、定位細胞胞器
- ✓ 可應用於細胞增生、存活、凋亡、週期...等檢測
- ✓ 可搭配螢光抗體進行多重染色 (multiplex)



細胞胞器、結構染劑

Product	Ex/Em (nm)	Cat. No.
Cell Plasma Membrane Staining Kit - Deep Red Fluorescence - Cytopainter	■ 640/660	ab219942
Cell Plasma Membrane Staining Kit - Orange Fluorescence - Cytopainter	■ 540/590	ab219941
ER Staining Kit - Green Fluorescence - Cytopainter	■ 441/551	ab139481
ER Staining Kit - Red Fluorescence - Cytopainter	■ 580/677	ab139482
Golgi Staining Kit - Green Fluorescence - Cytopainter	■ 473/534	ab139483
Golgi/ER Staining Kit - Cytopainter	■ Golgi 473/534 ■ ER 580/677	ab139485
Lysosomal Staining Kit - Blue Fluorescence - Cytopainter	■ 350/440	ab112135
Lysosomal Staining Kit - Deep Red - Cytopainter	■ 590/620	ab138896
Lysosomal Staining Kit - Green Fluorescence - Cytopainter	■ 490/525	ab112136
Lysosomal Staining Kit - Green Fluorescence (Ex405nm) - Cytopainter	■ 405/480	ab112138
Lysosomal Staining Kit - NIR Fluorescence - Cytopainter	■ 630/650	ab176746
Lysosomal Staining Kit - Orange - Cytopainter	■ 540/560	ab138895
Lysosomal Staining Kit - Red Fluorescence - Cytopainter	■ 575/600	ab112137
Lysosomal Staining Reagent - Blue - Cytopainter	■ 353/442	ab176825
Lysosomal Staining Reagent - Deep Red - Cytopainter	■ 596/619	ab176829
Lysosomal Staining Reagent - Green - Cytopainter	■ 450/505	ab176826
Lysosomal Staining Reagent - NIR - Cytopainter	■ 636/650	ab176824
Lysosomal Staining Reagent - Orange - Cytopainter	■ 542/556	ab176827
Lysosomal Staining Reagent - Red - Cytopainter	■ 575/597	ab176828
Lysosome/ER/Nuclear Staining Reagent - Cytopainter	■ Lysosome 568/667 ■ ER 440/565 ■ Nuclear 350/461	ab139488
Lysosome/Mitochondria/Nuclear Staining Reagent - Cytopainter	■ Lysosome 568/667 ■ Mitochondria 460/560 ■ Nuclear 350/491	ab139487
Mitochondrial Staining Kit - Blue - Cytopainter	■ 360/445	ab219939
Mitochondrial Staining Kit - Green Fluorescence - Cytopainter	■ 490/520	ab112143
Mitochondrial Staining Kit - NIR Fluorescence - Cytopainter	■ 640/660	ab176747
Mitochondrial Staining Kit - Orange (Ex405nm) - Cytopainter	■ 405/550	ab138897
Mitochondrial Staining Kit - Orange Fluorescence - Cytopainter	■ 545/575	ab112144
Mitochondrial Staining Kit - Red Fluorescence - Cytopainter	■ 580/600	ab112145

* 更多細胞螢光染劑產品請見內頁與背頁說明...





細胞結構染色劑

Product	Ex/Em (nm)	Cat. No.
Mitochondrial Staining Reagent - Blue - Cytopainter	■ 350/490	ab219940
Mitochondrial Staining Reagent - Green - Cytopainter	■ 498/520	ab176830
Mitochondrial Staining Reagent - NIR - Cytopainter	■ 640/659	ab176833
Mitochondrial Staining Reagent - Orange - Cytopainter	■ 545/575	ab176831
Mitochondrial Staining Reagent - Orange 405 - CytoPainter	■ 399/550	ab176834
Mitochondrial Staining Reagent - Red - Cytopainter	■ 575/600	ab176832
Nucleolar Staining Kit - Green Fluorescence - Cytopainter	■ 450/481	ab139475
Nucleolar/Nuclear Staining Kit - Green/ Red Fluorescence - Cytopainter	■ Nucleolus 450/490 ■ Nucleus 566/650	ab139474



纖維形肌動蛋白染色劑

Product	Ex/Em (nm)	Cat. No.
F-actin Staining Kit - Blue Fluorescence - Cytopainter	■ 350/450	ab112124
F-actin Staining Kit - Green Fluorescence - Cytopainter	■ 500/520	ab112125
F-actin Staining Kit - Orange Fluorescence - Cytopainter	■ 550/575	ab112126
F-actin Staining Kit - Red Fluorescence - Cytopainter	■ 594/610	ab112127
Phalloidin FITC Reagent	■ 492/518	ab235137
Phalloidin-iFluor 350 Reagent	■ 353/442	ab176751
Phalloidin-iFluor 405 Reagent	■ 400/421	ab176752
Phalloidin-iFluor 488 Reagent	■ 493/517	ab176753
Phalloidin-iFluor 514 Reagent	■ 520/547	ab176754
Phalloidin-iFluor 532 Reagent	■ 542/558	ab176755
Phalloidin-iFluor 555 Reagent	■ 556/574	ab176756
Phalloidin-iFluor 594 Reagent	■ 590/618	ab176757
Phalloidin-iFluor 633 Reagent	■ 634/649	ab176758
Phalloidin-iFluor 647 Reagent	■ 650/665	ab176759
Phalloidin-iFluor 680 Reagent	■ 683/700	ab176760
Phalloidin-iFluor 700 Reagent	■ 692/708	ab176761
Phalloidin-iFluor 750 Reagent	■ 752/778	ab176762
Phalloidin-iFluor 790 Reagent	■ 787/808	ab176763
Rhodamine Phalloidin Reagent	■ 546/575	ab235138

* 此系列產品僅適用於 Fixed cells



DNA 染色劑

Product	Ex/Em (nm)	Cat. No.
DAPI Staining Solution	■ DNA 358/461 ■ RNA 358/500	ab228549
DRAQ5™	■ 646/697	ab108410
DRAQ7™	■ 599/697	ab109202
Hoechst 33258 Staining Dye Solution	■ 352/461	ab228550
Hoechst 33342 Staining Dye Solution	■ 350/461	ab228551
Nuclear yellow (Hoechst S769121)	■ 355/495	ab138903
Propidium Iodide	■ 536/617	ab14083



細胞活性染劑

Product	Ex/Em (nm)	Cat. No.
Cell Viability Assay Kit (Fluorometric - Blue Ex 405 nm)	■ 405/460	ab176748
Cell Viability Assay Kit (Fluorometric - Blue)	■ 350/461	ab112120
Cell Viability Assay Kit (Fluorometric - Dual Green/Red)	■ Live cells 490/525 ■ Dead cells 540/620	ab112121
Cell Viability Assay Kit (Fluorometric - Green)	■ 490/520	ab112122
Fixable Cell Viability Assay Kit (Fluorometric - Blue Ex 405 nm) - Cytopainter	■ 410/450	ab176738
Fixable Cell Viability Assay Kit (Fluorometric - Blue) - Cytopainter	■ 353/442	ab176741
Fixable Cell Viability Assay Kit (Fluorometric - Deep Red) - Cytopainter	■ 649/660	ab176745
Fixable Cell Viability Assay Kit (Fluorometric - Green Ex 405 nm) - Cytopainter	■ 408/512	ab176739
Fixable Cell Viability Assay Kit (Fluorometric - Green) - Cytopainter	■ 498/521	ab176742
Fixable Cell Viability Assay Kit (Fluorometric - Orange Ex 405 nm) - Cytopainter	■ 398/550	ab176740
Fixable Cell Viability Assay Kit (Fluorometric - Orange) - Cytopainter	■ 547/573	ab176743
Fixable Cell Viability Assay Kit (Fluorometric - Red) - Cytopainter	■ 583/603	ab176744
Live and Dead Cell Assay	■ Live cells 494/515 ■ Dead cells 528/617	ab115347
Live Cell Labeling Kit - Blue Fluorescence - Cytopainter	■ 360/445	ab187963
Live Cell Labeling Kit - Blue Fluorescence (Ex405nm) - Cytopainter	■ 410/450	ab187966
Live Cell Labeling Kit - Green Fluorescence - Cytopainter	■ 495/515	ab187964
Live Cell Labeling Kit - Green Fluorescence (Ex405nm) - Cytopainter	■ 410/500	ab187967
Live Cell Labeling Kit - Red Fluorescence - Cytopainter	■ 649/659	ab187965



細胞週期、染色染劑

Product	Ex/Em (nm)	Cat. No.
Cell Cycle Assay Kit (Fluorometric - Green)	■ 490/520	ab112116
Cell Cycle Assay Kit (Fluorometric - Red)	■ 490/620	ab112117
Cell Proliferation Staining Reagent - Blue Fluorescence - Cytopainter	■ 403/454	ab176726
Cell Proliferation Staining Reagent - Deep Red Fluorescence - Cytopainter	■ 628/643	ab176736
Cell Proliferation Staining Reagent - Green Fluorescence - Cytopainter	■ 511/525	ab176735
Cell Proliferation Staining Reagent - Orange Fluorescence - Cytopainter	■ 542/556	ab176737
CFSE - Cell Labeling Kit	■ 492/517	ab113853
EdU Assay / EdU Staining Proliferation Kit (iFluor 488)	■ 491/520	ab219801
EdU Assay / EdU Staining Proliferation Kit (iFluor 647)	■ 649/664	ab222421
Propidium Iodide Flow Cytometry Kit	■ 536/617	ab139418



細胞追蹤染劑

Product	Ex/Em (nm)	Cat. No.
Cell Tracking Dye Kit - Blue - Cytopainter	■ 360/445	ab138890
Cell Tracking Dye Kit - Deep Red - Cytopainter	■ 630/650	ab138894
Cell Tracking Dye Kit - Green - Cytopainter	■ 490/520	ab138891
Cell Tracking Dye Kit - Orange - Cytopainter	■ 540/560	ab138892
Cell Tracking Dye Kit - Red - Cytopainter	■ 575/600	ab138893
Cell Tracking Red Dye Kit - Longer cell staining, DMSO-free	■ 554/575	ab269446



粒檢 線測 體試 膜劑 電位

Product	Ex/Em (nm)	Cat. No.
JC-1 - Mitochondrial Membrane Potential Assay Kit	■ Monomer 475/530 ■ Aggregate 475, 535/590	ab113850
JC-10 Mitochondrial Membrane Potential Assay Kit (Flow Cytometry)	■ Monomer 490/520 ■ Aggregate 490/570	ab112133
JC-10 Mitochondrial Membrane Potential Assay Kit (Microplate)	■ Monomer 490/520 ■ Aggregate 490/570	ab112134
NIR Mitochondrial Membrane Potential Assay Kit (Microplate)	■ 640/680	ab112150
Orange Mitochondrial Membrane Potential Assay Kit (Flow Cytometry)	■ 540/590	ab138898
Orange Mitochondrial Membrane Potential Assay Kit (Microplate)	■ 540/590	ab138899
TMRE-Mitochondrial Membrane Potential Assay Kit	■ 549/575	ab113852
TMRM Assay Kit (Mitochondrial Membrane Potential)	■ 548/573	ab228569



鈣離子 流動 檢測 試劑

Product	Ex/Em (nm)	Cat. No.
Calcium Flux Assay Kit (Flow cytometry)	■ 490/525	ab233472
Fluo-2 AM, green fluorescent Ca ²⁺ binding dye	■ 490/520	ab142775
Fluo-3 AM, Fluorescent calcium Indicator	■ 506/526	ab145254
Fluo-4 AM, Fluorescent labeling reagent	■ 494/516	ab241082
Fluo-4 Assay Kit (Calcium)	■ 490/525	ab228555
Fluo-5 AM, Fluorescent labeling reagent (Fluo-4 analogue)	■ 494/516	ab241083
Fluo-8 AM, green fluorescent calcium binding dye	■ 490/520	ab142773
Fluo-8 Calcium Assay Kit - Medium Removal	■ 490/525	ab112128
Fluo-8 Calcium Flux Assay Kit - No Wash	■ 490/525	ab112129
Fura-2 AM, Ca ²⁺ selective fluorescent indicator	■ 340, 380/505	ab120873
Fura-2 Calcium Flux Assay Kit - No Wash, Ratiometric	■ 340, 380/510	ab176766
Indo-1 AM, fluorescent Ca ²⁺ indicator	■ 346/475, 405	ab142778
Rhod-2 AM, fluorescent Ca ²⁺ indicator	■ 552/581	ab142780
Rhod-4 Calcium Assay kit - Medium Removal No Wash	■ 540/590	ab112156
Rhod-4 Calcium Assay Kit - No Wash	■ 540/590	ab112157



神影 經像 科探 學針

Product	Ex/Em (nm)	Cat. No.
CRANAD-2, Near-infrared Amyloid-beta fluorescent probe	■ 640/715	ab141775
MCAAD-3, Near-infrared probe for imaging of amyloid plaques	■ 596/685	ab216983



封片膠

Product	Cat. No.
Anti-Fade Fluorescence Mounting Medium - Aqueous, Fluoroshield	ab104135
Glycerol Mounting Medium - Anti-Fade with DAPI	ab188804
Mounting Medium With DAPI - Aqueous, Fluoroshield	ab104139
Mounting Medium With PI - Aqueous, Fluoroshield	ab104129

產品資訊 · 促銷優惠

請洽 Abcam 台灣代理 [伯森生技](#)



伯森生技

伯森生技客服專線

北部 0800-059668 南部 0800-259988

